

## Adapt Lubricating Deodorant

### General Information

- **Adapt** lubricating deodorant neutralizes odors, eases emptying, and helps prevent pouch static and sticking
- It is not a perfume or masking agent
- Adapt lubricating deodorant should be used in drainable and closed pouching systems only

### For Drainable Pouches

1. Apply pouch.
2. Bottle: add approximately 5 mL or 1 tsp to bottom of pouch (Figure 1).  
Packet: add the complete contents of the packet to bottom of pouch (Figure 2). Close the pouch in the normal manner.
3. Rub to thoroughly coat the inside of the pouch (Figure 3). Repeat Step 2 every time you empty your pouch.

### For Closed Pouches

1. Add Adapt lubricating deodorant before removing the release liner on your skin barrier.
2. Bottle: add approximately 5 mL or 1 tsp through pouch opening (Figure 4).  
Packet: add the complete contents of the packet through pouch opening.
3. Apply pouch.
4. Rub to thoroughly coat the inside of the pouch (Figure 3).



### FREQUENTLY ASKED QUESTIONS

Does it matter if the Adapt lubricating deodorant comes in contact with the filter on the pouch?

When coating the inside of the pouch, it is best to avoid the filter.

Is Adapt lubricating deodorant covered by Medicare?\*

Yes, Adapt lubricating deodorant is billed per ounce using HCPCS code A4394 (Ostomy deodorant, with or without lubricant, for use in ostomy pouch, liquid). Medicare has not published a usual maximum quantity for this code. Both the 8 oz bottle and 8 mL packets (50 packets = 13.5 oz, 1 box = 13.5 oz) are billed using the same code. Check with your supplier to confirm.

Do I need to use the Adapt lubricating deodorant every time I empty my pouch?

Yes. For this ostomy accessory to work effectively, it must be added each time your pouch is drained.

Can I just put a few drops in the bottom of my pouch and not rub it against the sides of the pouch?

No, this is not recommended since it will not be effective as a pouch lubricant or as a pouch deodorant. Rubbing the pouch helps to coat the inside of the pouch film.

Isn't it easier to use a cooking spray or oil instead to coat the inside of my pouch?

This is not recommended because cooking oil may cause the seal of your pouch to fail. Also, cooking oil sprays do not eliminate odor.

## Adapt Lubricating Deodorant

### FREQUENTLY ASKED QUESTIONS

<p>The 8 oz bottle is not convenient to take into a public restroom. Is there a smaller bottle available?</p>	<p>Adapt lubricating deodorant is also available in single-use packets. This is much more discrete and convenient when you are away from home.</p>
<p>It doesn't have a scent. How do I know it is working?</p>	<p>Adapt lubricating deodorant is not a perfume or a fragrance used to mask or cover-up odors. It does not have any odor by itself, but it will neutralize odors in the pouch.</p>
<p>Can the Adapt lubricating deodorant come in contact with the stoma?</p>	<p>Yes, it has been tested and it will not harm your stoma.</p>
<p>What if Adapt lubricating deodorant comes in contact with my skin?</p>	<p>This product is not a skin irritant. It is water-soluble, so it will easily be removed from the skin with a moist cloth. Be sure to wash it off your skin before applying your pouch, so it does not interfere with adhesion.</p>



#### Adapt lubricating deodorant

Stock No 78500  
(8 oz bottle)

Stock No 78501  
(box of 50, 8 mL packets)

For product questions, sampling needs, or detailed clinical questions concerning our products in the US, call **1.888.808.7456**. In Canada call **1.800.263.7400**.

Routine follow-up with your healthcare professional is recommended.

Prior to using any ostomy products/accessories be sure to read all product inserts and labels.

\*The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, and can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carrier or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

The Hollister logo, Adapt and "Healthy skin. Positive outcomes." are trademarks of Hollister Incorporated.  
© 2017 Hollister Incorporated. 923061-1217 US-00214

**Hollister Incorporated**  
2000 Hollister Drive  
Libertyville, Illinois 60048  
1.800.323.4060

[www.hollister.com](http://www.hollister.com)

**Hollister Limited**  
95 Mary Street  
Aurora, Ontario L4G 1G3  
1.800.263.7400